


Elements	Presenter	Primary experimental data
Transcripts (mRNAs, ncRNAs, TSS, UTRs, miRNAs)	Celniker Lai	Tiling arrays, RNASeq, RT-PCR/RACE, Mass Spec, 3'UTR clone library, UAS-miRNA flies.
Enhancers, promoters, silencers, insulators	White	ChIP/chip, ChIP/seq, TF-tagged strains, anti-TF antibodies.
Chromatin marks	Pirrotta Ahmad	ChIP/chip and ChIP/seq of chromosome-associated proteins and nucleosome
DNA Replication origins	Orr-Weaver	ChIP/chip and ChIP/seq of essential initiator proteins, origin mapping and DNA copy number in differentiated tissues
Data Coordination Center	Suzi Lewis	www.modencode.org

A cis-regulatory map of the *Drosophila* genome


- **Data generation**
 - Kevin White, U. Chicago (Antibody pipeline, ChIP chip pipeline)
 - Bing Ren, UCSD (Ab validation, ChIP chip pipeline)
- **Computational identification of cis-regulatory sequences**
 - Manolis Kellis, MIT (Comparative motif analysis, ChIP chip analysis)
- **Validation**
 - Jim Posakony, UCSD (Promoters/Enhancers)
 - Steve Russell, U. Cambridge (Insulators/Silencers)
 - Hugo Bellen, Baylor (element “necessity” validations)

Cis regulatory Element types


– Promoters




– Enhancers



– Silencers and Insulators



A Cis-Regulatory Annotation Pipeline



Target Proteins to define each type of Cis regulatory sequence

– Promoters

- **Basal transcriptional machinery (e.g. Pol II)**
- **Histone state (e.g. H3K4 tri-methylation)**

– Enhancers


- **Site-specific TFs (target >350)**
- **Histone state (e.g. H3K4 mono-methylation)**
- **Co-factors (e.g. p300)**

– Silencers and Insulators

- **Silencers – (e.g. PcG)**
- **Insulators - (e.g. Su(Hw), Mod(mdg4))**





What samples to ChIP?


Cumulative Gene Expression During Development (by molecular classification)



Data from Arbeitman et al. Science 2002

Preliminary data - Chromatin signatures (B Ren)





Drosophila S2 Clusters

❖ Active Promoters (3997)

- + H3K4me3
- + H3K9Ac
- + H3K18Ac
- + H3K27Ac
- + RNAP
- H3
- H3K4me1

❖ Enhancers (4242)

- + H3K4me1*
- + H3K18Ac
- + H3K27Ac
- H3
- H3K4me3
- H3K9Ac
- RNAP


Sample Collections Aug-Nov 2007

<i>Stages</i>	<i>Amt. of Bio. Mat. (g)</i>	<i>Chromatin Equivalent</i>
0-4h	0.994	3300ug
4-8h	0.159	530ug
12-16h	0.762	2540ug
16-20h	0.681	2270ug
20-24h	1.077	3590ug
L1	1.098	3660ug
L2	0.880	2930ug
L3	1.890	6300ug
Pupae	1.090	3630ug
Female Adults	0.757	2520ug
Male Adults	0.460	1530ug

Histone Modifications
H3K4me1
H3K4me3
H3K9Ac
H3K9me3
H3K27Ac
H3K27me3
(PolII)

RNA Extractions (for the moment homogenized in TRIzol)

<i>Date</i>	<i>Sample</i>	<i>Amt of Bio. Material</i>
7.25.07	Pupae	110.8mg
9.10.07	L3	135mg/98mg
9.20.07	0-4h	195mg/94mg
9.21.07	4-8H	329mg
10.10.07	20-24h	55mg
10.11.07	16-20h	110mg




Top: ChIP using 100ug adult female chromatin (30mg Bio. Material) and 10ug H3K4Me3 antibody, DNA was amplified by LM-PCR and labeled with the Nimblegen protocol (end labeled random priming).

Middle: Same as above except using 300ug adult female chromatin.

Bottom: ChIP using 100ug adult female chromatin (30mg Bio. Material) and 10ug H3K4Me3 antibody, DNA was amplified by WGA2 and labeled with random priming.


ChIP milestones

- **year 1 – Genome-wide scans with existing Abs to basal TF machinery and chromatin marks; 300 experiments total**
 - Chromatin time course plus two Cell lines (S2, Kc) = 909 arrays (progress: 315 complete; 72-96 per week)
 - **year 2 – 1,000 total arrays (+700)**
 - **year 3 – 2,000 total arrays (+1,000)**
 - **year 4 – 3,500 total arrays (+1,500)**
- 
- Site-specific TFs**

Technical validations:

- **Replicate experiments**
- **Independent Antibodies (when available)**
- **qPCR/Solexa**


Regulatory motif discovery pipeline (M. Kellis)




Validation milestones (J. Posakony)

- **year 1 – Establish criteria for choosing elements to validate, test and characterize in vivo reporter gene system, and train new personnel**
- **year 2 - Test 100 total elements each of enhancers and promoters (Posakony), 20 each of silencers and insulators (Russell)**
- **year 3 – Test 200 total elements each of enhancers and promoters, 40 each of silencers and insulators**
- **year 4 – Test 300 total elements each of enhancers and promoters, 60 each of silencers and insulators**

Reporter Construct Design



A close-up photograph of two young, light brown rabbits with long, upright ears. They are nestled among green grass and numerous white daisies with yellow centers. The rabbit on the left is looking slightly towards the camera, while the one on the right is looking down. The background is a soft-focus field of more daisies.





Antibody production pipeline

- **Protein domains (~100a.a.) 35 antisera for 22 factors yielded 15 validated antibodies (43%) for 13 factors (59%)**
- **100 peptide based antibodies for 50 factors; 81 antibodies tested by microwestern (w/Rich Jones and Mark Ciaccio) 18/81 show specific bands in WB = 22%**
- **104 protein domains injected in December; 100 peptides representing 50 factors injected at beginning of Feb**
- **100+ domains have completed production phase, awaiting injection**

Antibody Production Milestones

Expectation that 20% of antibodies will be useful for ChIP

- **year 1 – produce 350 total antibodies**
 - **Progress (funding started 8/07): 339 injected with 100 more protein domains this month, 116 antisera received.**
Validated *in situ*/Western: 33/116 = 28%
- **year 2 – 850 total antibodies produced (+500)**
- **year 3 – 1350 total antibodies produced (+500)**
- **year 4 – 1750 total antibodies produced (+400)**

Total proteins ~500 (2-3 domains each)

**Total proteins that are ChIP validated ~350
(40-70% with multiple ChIP validated Abs)**

Please help us prioritize Ab and ChIP experiments

- **Vote for your favorite
transcription factor!!!**
[**http://www.modencode.org/Vot**](http://www.modencode.org/Vot)
- **Send us your antibodies to test,
validate and post ChIP data
online!!! Nicolas Negre**
[**nnegre@bsd.uchicago.edu**](mailto:nnegre@bsd.uchicago.edu)